

# **CAPITULOVIII**

## **ANEXOS**

**ANEXO 1. PROSPECTO PARA LA DETERMINACION DE SARS-COV-2 (COVID-19)**



**NovaLisa®**

**SARS-CoV-2 (COVID-19) IgM**

**ELISA**

**For Performance Study ONLY**

English.....	2
Deutsch.....	7
Bibliography / Literatur / Bibliographie / Bibliografía / Bibliografia / Bibliografía.....	14
Abbreviations / Abkürzungen / Abréviations / Abbreviazioni / Abreviaciones / Abreviações.....	14
Symbols Key / Symbolschlüssel / Explication des Symboles / Legenda / Símbolos / Tabela de símbolos.....	15
Summary of Test Procedure / Kurzanleitung Testdurchführung / Résumé de la procédure de test / Schema della procedura / Resumen de la técnica / Resumo do Procedimento de Teste.....	15

---

**Product Number: COVM0940 (96 Determinations)**

---

## ENGLISH

### 1. INTRODUCTION

End of 2019, a novel respiratory disease emerged in the city of Wuhan, Hubei Province of the People's Republic of China, and soon spread rapidly within the country and worldwide. The causative agent was identified as severe acute respiratory syndrome coronavirus 2 (SARS-CoV-2). SARS-CoV-2 (2019-nCoV), like the closely related SARS coronavirus (SARS-CoV), belongs to the genus Betacoronavirus within the family of coronaviruses. The zoonotic reservoir of the virus appears to be bats.

Coronaviruses are enveloped, positive single-stranded large RNA viruses that infect humans, but also a wide range of animals. The common human coronaviruses NL63, 229E, OC43 and HKU1 are widespread especially throughout the winter months. They are responsible for up to one third of all acute respiratory diseases, typically with mild symptoms (common cold). More than 80 % of the adult population have antibodies against human coronaviruses. The immunity from previous infections lasts only for a short period of time. Therefore, reinfections with the same pathogen are possible just after one year.

SARS-CoV-2 is predominantly transmitted by droplet infection via coughing or sneezing and through close contact with infected patients. In theory, smear infection and infection through the conjunctiva of the eyes are also possible.

The incubation period is in the median 5-6 days (and up to 14 days maximum).

The clinical manifestations of SARS-CoV-2-related COVID-19 disease include fever, cough, respiratory problems and fatigue. In most patients the infection manifests with symptoms of a mild febrile illness with irregular lung infiltrates.

The initial clinical sign of COVID-19 which allowed case detection was pneumonia. But it turned out that the course of the disease is non-specific and varies widely, from asymptomatic courses to severe pneumonia with lung failure and death. However, based on current knowledge, around 80 % of the illnesses are mild to moderate.

Although severe courses of the disease also occur in younger patients and people without previous illness, the following groups of people have an increased risk of serious forms of the disease: elderly people (with a steadily increasing risk from around 50-60 years of age), smokers and people with certain diseases of the cardiovascular system or the lungs, patients with chronic liver diseases, diabetes mellitus, cancer, or patients with a weakened immune system (e.g. due to immune deficiencies or by taking drugs that suppress the immune system).

Currently, there is no specific treatment or vaccine available against SARS-CoV-2 infection.

Species	Disease	Symptoms e.g.	Transmission route
SARS-CoV-2 (severe acute respiratory syndrome coronavirus 2)	COVID-19	the course of the disease is unspecific, diverse and varies greatly, from asymptomatic courses to severe pneumonia with lung failure and death	primary mode of transmission: droplet infection; smear infections and infections via the conjunctiva of the eyes are theoretically possible

The presence of pathogen or infection may be identified by

- Nucleic acid testing (NAT): e.g. RT-PCR
- Serology: detection of antibodies by ELISA, immunoblot

### 2. INTENDED USE

The SARS-CoV-2 (COVID-19) IgM ELISA is intended for the qualitative determination of IgM class antibodies against SARS-CoV-2 in human serum or plasma (citrate, heparin) to support the diagnosis of COVID-19 disease and constitutes a supplement to direct pathogen detection.

### 3. PRINCIPLE OF THE ASSAY

The qualitative immunoenzymatic determination of specific antibodies is based on the ELISA (Enzyme-linked Immunosorbent Assay) technique.

Microtiterplates are coated with specific antigens to bind corresponding antibodies of the sample. After washing the wells to remove all unbound sample material a horseradish peroxidase (HRP) labelled conjugate is added. This conjugate binds to the captured antibodies. In a second washing step unbound conjugate is removed. The immune complex formed by the bound conjugate is visualized by adding Tetramethylbenzidine (TMB) substrate which gives a blue reaction product.

The intensity of this product is proportional to the amount of specific antibodies in the sample. Sulphuric acid is added to stop the reaction. This produces a yellow endpoint colour. Absorbance at 450/620 nm is read using an ELISA Microtiterplate reader.

## 4. MATERIALS

### 4.1. Reagents supplied

- Microtiterplate: 12 break-apart 8-well strip-off strips coated with SARS-CoV-2 antigens, in resealable aluminium foil.
- IgM Sample Dilution Buffer: 1 bottle containing 100 mL of phosphate buffer (10 mM) for sample dilution; pH 7.2 ± 0.2; anti-human IgG (RF Absorbent) coloured green; ready to use; white cap; ≤ 0.0015 % (w/v) CHSE/MT (3.1).
- Stop Solution: 1 bottle containing 15 mL sulphuric acid, 0.2 mol/L; ready to use; red cap.
- Washing Buffer (20x conc.): 1 bottle containing 50 mL of a 20-fold concentrated phosphate buffer (0.2 M) pH 7.2 ± 0.2, for washing the wells; white cap.
- Conjugate: 1 bottle containing 20 mL of peroxidase labelled antibody to human IgM in phosphate buffer (10 mM) coloured red; ready to use; black cap.
- TMB Substrate Solution: 1 bottle containing 15 mL 3,3',4,5-tetramethylbenzidine (TMB), = 0.1 %; ready to use; yellow cap.
- Positive Control: 1 vial containing 2 mL control, coloured yellow; ready to use; red cap; ≤ 0.02 % (w/v) MT.
- Cut-off Control: 1 vial containing 2 mL control, coloured yellow; ready to use; green cap; ≤ 0.02 % (w/v) MT.
- Negative Control: 1 vial containing 2 mL control, coloured yellow; ready to use; blue cap; ≤ 0.0015 % (w/v) CHSE/MT (3.1).

For hazard and precautionary statements see 12.1.

For potential hazardous substances please check the safety data sheet.

### 4.2. Materials supplied

- 1 Cover foil
- 1 Instruction for use (IFU)
- 1 Plate layout

### 4.3. Materials and Equipment needed

- ELISA Microtiterplate reader, equipped for the measurement of absorbance at 450/620 nm.
- Incubator 37 °C
- Manual or automatic equipment for reading Microtiterplate wells
- Pipettes to deliver volumes between 10 and 1000 µL
- Vortex tube mixer
- Distilled water
- Disposable tubes

## 5. STABILITY AND STORAGE

Store the kit at 2...8 °C. The opened reagents are stable up to the expiry date stated on the label when stored at 2...8 °C.

## 6. REAGENT PREPARATION

It is very important to bring all reagents and samples to room temperature (20...25 °C) and mix them before starting the test run!

### 6.1. Microtiterplate

The break-apart strip-off strips are coated with SARS-CoV-2 antigens. Immediately after removal of the strips, the remaining strips should be resealed in the aluminium foil along with the desiccant supplied and stored at 2...8 °C.

### 6.2. Washing Buffer (20x conc.)

Dilute Washing Buffer 1 × 10, e.g. 10 mL Washing Buffer + 190 mL distilled water. The diluted buffer is stable for 5 days at room temperature (20...25 °C). In case crystals appear in the concentrate, warm up the solution to 37 °C e.g. in a water bath. Mix well before dilution.

### 6.3. TMB Substrate Solution

The reagent is ready to use and has to be stored at 2...8 °C, away from the light. The solution should be colourless or could have a slight blue tinge. If the substrate turns into blue, it may have become contaminated and should be thrown away.

## 7. SAMPLE COLLECTION AND PREPARATION

Use human serum or plasma (oblate, heparin) samples with this assay. If the assay is performed within 5 days after sample collection, the samples should be kept at 2...8 °C; otherwise they should be aliquoted and stored deep frozen (-70...-20 °C). If samples are stored frozen, mix thawed samples well before testing. Avoid repeated freezing and thawing. Heat inactivation of samples is not recommended.

### 7.1. Sample Dilution

Before assaying, all samples should be diluted 1 × 100 with IgM Sample Dilution Buffer. Dispense 10 µL sample and 1 mL IgM Sample Dilution Buffer into tubes to obtain a 1:100 dilution and thoroughly mix with a Vortex.

---

## 8. ASSAY PROCEDURE

Please read the instruction for use carefully before performing the assay. Result reliability depends on strict adherence to the instruction for use as described. The following test procedure is only validated for manual procedure. If performing the test on ELISA automatic systems we recommend increasing the washing steps from three up to five and the volume of Washing Buffer from 300 µL to 350 µL to avoid washing effects. Pay attention to chapter 12. Prior to commencing the assay, the distribution and identification plan for all samples and standards/controls (duplicates recommended) should be carefully established on the plate layout supplied in the kit. Select the required number of microtiter strips or wells and insert them into the holder.

Perform all assay steps in the order given and without any delays.

A clean, disposable tip should be used for dispensing each standard/control and sample.

Adjust the incubator to  $37 \pm 1$  °C.

1. Dispense 100 µL standards/controls and diluted samples into their respective wells. Leave well A1 for the Substrate Blank.
2. Cover wells with the foil supplied in the kit.
3. Incubate for 1 hour  $\pm$  5 min at  $37 \pm 1$  °C.
4. When incubation has been completed, remove the foil, aspirate the content of the wells and wash each well three times with 300 µL of Washing Buffer. Avoid overflows from the reaction wells. The interval between washing and aspiration should be > 5 sec. At the end carefully remove remaining fluid by tapping strips on tissue paper prior to the next step.  
Note: Washing is important! Insufficient washing results in poor precision and false results.
5. Dispense 100 µL Conjugate into all wells except for the Substrate Blank well A1.
6. Incubate for 30 min at room temperature (20...25 °C). Do not expose to direct sunlight.
7. Repeat step 4.
8. Dispense 100 µL TMB Substrate Solution into all wells.
9. Incubate for exactly 15 min at room temperature (20...25 °C) in the dark. A blue colour occurs due to an enzymatic reaction.
10. Dispense 100 µL Stop Solution into all wells in the same order and at the same rate as for the TMB Substrate Solution, thereby a colour change from blue to yellow occurs.
11. Measure the absorbance at 450/620 nm within 30 min after addition of the Stop Solution.

### 8.1. Measurement

Adjust the ELISA Microtiterplate reader to zero using the Substrate Blank.

If - due to technical reasons - the ELISA Microtiterplate reader cannot be adjusted to zero using the Substrate Blank, subtract its absorbance value from all other absorbance values measured in order to obtain reliable results!

Measure the absorbance of all wells at 450 nm and record the absorbance values for each standard/control and sample in the plate layout.

Dichromatic measurement using a reference wavelength of 620 nm is recommended.

Where applicable calculate the mean absorbance values of all duplicates.

---

## 9. RESULTS

### 9.1. Run Validation Criteria

In order for an assay to be considered valid, the following criteria must be met:

- **Substrate Blank:** Absorbance value < 0.100
- **Negative Control:** Absorbance value < 0.200 and < Cut-off
- **Cut-off Control:** Absorbance value 0.150 – 1.300
- **Positive Control:** Absorbance value > Cut-off

If these criteria are not met, the test is not valid and must be repeated.

### 9.2. Calculation of Results

The Cut-off is the mean absorbance value of the Cut-off Control determinations.

Example: Absorbance value Cut-off Control 0.44 + absorbance value Cut-off control 0.42 =  $0.86 / 2 = 0.43$   
Cut-off = 0.43

#### 9.2.1. Results in Units [NTU]

Sample (mean) absorbance value  $\times 10$  = [NovaTec Units = NTU]  
Cut-off

Example:  $\frac{1.591 \times 10}{0.43} = 37$  NTU (Units)

---

### 9.3. Interpretation of Results

Cut-off	ID NTU	
Positive	> 11 NTU	Antibodies against the pathogen are present. There has been a contact with the antigen (pathogen resp. vaccine).
Equivocal	9 – 11 NTU	Antibodies against the pathogen could not be detected clearly. It is recommended to repeat the test with a fresh sample in 2 to 4 weeks. If the result is equivocal again the sample is judged as negative.
Negative	< 9 NTU	The sample contains no antibodies against the pathogen. A previous contact with the antigen (pathogen resp. vaccine) is unlikely.

Diagnosis of an infectious disease should not be established on the basis of a single test result. A precise diagnosis should take into consideration clinical history, symptomatology as well as serological data.  
In immunocompromised patients and newborns serological data only have restricted value.

#### 9.3.1. Antibody isotypes and State of Infection

Serology	Significance
IgM	Characteristic of the primary antibody response High IgM titer: => suggests a current or very recent infection
IgG	Follows IgM production Characteristic of the secondary antibody response May persist for several years High IgG titer with low IgM titer: => may indicate past infection
IgA	Produced in mucosal linings throughout the body (=> protective barrier) Usually produced early in the course of the infection

### 10. SPECIFIC PERFORMANCE CHARACTERISTICS

The results refer to the groups of samples investigated; these are not guaranteed specifications.

For further information about the specific performance characteristics please contact NovaTec Immundiagnostica GmbH.

#### 10.1. Precision

Intraassay	n	Mean (E)	CV (%)
#1	24		
#2	24		
#3	24		

Interassay	n	Mean (NTU)	CV (%)
#1	12		
#2	12		
#3	12		

#### 10.2. Diagnostic Specificity

The diagnostic specificity is defined as the probability of the assay of scoring negative in the absence of the specific analyte. It is xx.xx % (95% confidence interval: xx.xx % - xx.xx %).

#### 10.3. Diagnostic Sensitivity

The diagnostic sensitivity is defined as the probability of the assay of scoring positive in the presence of the specific analyte. It is xx.xx % (95% confidence interval: xx.xx % - xx.xx %).

#### 10.4. Interferences


#### 10.5. Cross Reactivity

### 11. LIMITATIONS OF THE PROCEDURE

Bacterial contamination or repeated freeze-thaw cycles of the sample may affect the absorbance values.


## ANEXO 2 ESTRUCTURA DE LOS CORONAVIRUS

Figura 1




Estructura genérica de los coronavirus (figura extraída de Wikipedia Commons)

Figura 2



# ANEXO 5 CICLO DE REPLICACION DEL SARS-COV-2

Figura 3




Ciclo de replicación del SARS-CoV-2 (Imagen modificada de Vega Asentio e Ignacio Lopez-Goñi, bajo licencia CC BY-SA).



## ANEXO 4 MECANISMO DE TRANSMISIÓN DEL SARS-COV-2

Figura 4. Emisión de secreciones respiratorias de diferentes tamaños por parte del caso índice (emisor) a un contacto expuesto (receptor) a una distancia menor a dos metros, con y sin mascarilla



Puntos azules: gotas balísticas ( $\geq 100 \mu\text{m}$ ); puntos verdes: aerosoles grandes ( $>15 \mu\text{m}$  hasta  $100 \mu\text{m}$ ); puntos naranjas: aerosoles intermedios ( $>5 \mu\text{m}$  hasta  $15 \mu\text{m}$ ); puntos rojos: aerosoles de pequeño tamaño ( $\leq 5 \mu\text{m}$ )

TOMA DE MUESTRA PARA PCR-RT

PARA SARS-COV2




# ANEA DE PROCESAMIENTO DE LAS MUESTRAS

## TECNICA DE ELISA PARA SARS-COV-2

### PARA SARS-COV2

1. Los micropocillos están recubiertos con anticuerpos específicos para la sustancia de interés.




Cuanto menos color azul, o más rojo, más sustancia de interés detectada.


## ANEXO 9

### TECNICA DE ELISA PARA SARS-COV-2


1. Los micropocillos están recubiertos con anticuerpos específicos para la sustancia de interés.




2. El conjugado compite con la sustancia/controles de interés para los sitios de unión de anticuerpos.



3. El conjugado y sustancia/controles de interés permanecen unidos en los pocillos.



4. Se añade sustrato para producir un cambio de color.



5. Los resultados se leen visualmente o en un lector – cuanto menos color azul, o más rojo, más sustancia de interés detectada.



## ANEXO 10.INTERPRETACIONDE RESULTADOS COMBINANDO PCR- RT Y DETECCION DE

PCR	Ag	IgM	IgG	Interpretación
+	-	-	-	Fase presintomática
+	+/-	+/-	+/-	Fase inicial (aprox. 1-7 días)
+/-	-	+	+/-	2ª Fase (8-14 días)
+/-	-	++	++	3ª Fase > 15 días
-	-	+/-	++	Infección pasada (inmune)

### ANTICUERPOS